# PNI WFO175

Oxygen facial beauty machine / Aparat de hidratare faciala

EN	User manual .....	3
BG	Ръководство за употреба .....	5
DE	Benutzerhandbuch .....	7
ES	Manual de usuario .....	9
FR	Manuel utilisateur .....	11
HU	Használati utasítás .....	13
IT	Manuale utente .....	15
NL	Handleiding .....	17
PL	Instrukcja obsługi .....	19
RO	Manual de utilizare .....	21

## Charging the battery

Before the first use, please charge the battery with the USB power cable included in the package. Do not use the oxygen facial beauty machine while it is charging.

## Usage

Fill the tank with still mineral water, distilled water or cosmetic water. Do not use concentrated skin cleansing lotions, dilute them with water. Fix the tank well on the device. Press the ON/OFF button to turn on the oxygen facial beauty machine. Press and hold the ON/OFF button to start spraying. Release the button to stop spraying. Briefly press the ON/OFF button to turn off the oxygen facial beauty machine.

This device is recommended for deep cleaning and skin hydration. Thus, you will have softer, brighter and more elastic skin.

## Safety warnings

- Clean the oxygen facial beauty machine with a slightly damp cloth, then wipe it well with a dry cloth.
- Keep the oxygen facial beauty machine in well-ventilated areas, away from excessive heat and humidity.
- In case of malfunction, do not try to repair the device yourself. Call a specialized service center.
- Do not leave the oxygen facial beauty machine within the reach

of children.

- Do not submerge the oxygen facial beauty machine in water.
- If you feel pain or discomfort, stop using this device immediately.
- Do not use this oxygen facial beauty machine if you have wounds or bruises on your face, if you are allergic to metals, if you suffer from skin diseases, if you are pregnant, if you suffer from lung diseases.
- If you suffer from serious diseases, consult your doctor before using this device.
- To charge the oxygen facial beauty machine, use the USB cable included in the package.
- Use only the accessories included in the package.

## Technical specifications

- Material: aluminum alloy, copper and plastic
- Power supply: 5 V / 0.5 A (via USB C port)
- Rated voltage: 7.4 V
- Battery capacity: 800mAh
- Rated current: 1.8A
- Power: max. 4.8W
- Tank volume: 20 ml / 40 ml
- Dimensions: 110 x 45 x 170 mm
- Operating temperature: -26°C ~ +80°C

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